



CCI Vials, Caps & Septas ...



◆ OEM based products. ◆ Certificate of analysis comes with all materials. ◆ Quality guaranteed. ◆ Internationally accepted

Chromachemie Laboratory Private Limited introduces "C.C.I" vials and closures. CCI vials/Caps are manufactured under stringent quality procedures to provide better quality products at an affordable price. Vials/Closures are HPLC grade and USP certified.

Do you know vials can add contamination to your sample?
Kindly Consult us for selecting right vials for your applications.

ND 9 - Screw Neck Short thread for HPLC. [1.5-2ml of 11.6mm x 32mm]

Vials



CCI - V917
1.5ml wide opening short screw-thread vial, clear



CCI - V927
1.5ml wide opening short screw-thread vial with write-on spot, clear



CCI - V935
1.5ml wide opening short screw-thread vial, amber



CCI - V945
1.5ml wide opening short screw-thread vial with write-on spot, amber




CCI - VP91
0.3ml PP Short Screw-Thread Micro-Vial, clear




CCI - VP93
0.3ml PP Short Screw-Thread Micro-Vial, amber


Caps / Septas




CCI - SC 9291
red PTFE/white Silicone septa, 9mm blue short screw-thread polypropylene cap, 6mm centre hole



CCI - SC 9491
blue PTFE/white Silicone septa, 9mm blue short screw-thread polypropylene cap, 6mm centre hole



CCI - SC 92291
pre-slit red PTFE/white Silicone septa, 9mm blue short screw-thread polypropylene cap, 6mm centre hole



CCI - SC 944491
Pre-slit blue PTFE/white Silicone septa, 9mm blue short screw-thread polypropylene cap, 6mm centre hole

ND 11 - Crimp type for HPLC / GC. [1.5-2ml of 11.6mm x 32mm]

Vials



CCI - V1117
1.5ml wide opening crimp-top vial, clear



CCI - V1127
1.5ml wide opening crimp-top vial with write-on spot, clear




CCI - V1135
1.5ml wide opening crimp-top vial, amber




CCI - V1145
1.5ml wide opening crimp-top vial with write-on spot, amber

Caps / Septas



CCI - SC 112111
red PTFE/white Silicone septa, 11mm crimp-top aluminum cap, 5.5mm centre hole



CCI - S 113
red PTFE/white Silicone/red PTFE septa



CCI - C 111
crimp-top aluminum cap, 5.5mm centre hole



CCI - SC 1122111
pre-slit red PTFE/white Silicone septa, 11mm crimp-top aluminum cap, 5.5mm centre hole

ND 11 - Snap type for HPLC. [1.5-2ml of 11.6mm x 32mm]

Vials



CCI - V1017
1.5ml wide opening snap-top vial, clear



CCI - V1027
1.5ml wide opening snap-top vial with write-on spot, clear



CCI - V1035
1.5ml wide opening snap-top vial, amber



CCI - V1045
1.5ml wide opening snap-top vial with write-on spot, amber




CCI - Vp101
0.3ml PP Snap Ring Micro-Vial, clear




CCI - Vp103
0.3ml PP Snap Ring Micro-Vial, amber


Caps / Septas




CCI - SC 102101
red PTFE/white Silicone septa, 11mm natural snap-top polypropylene cap, 6mm centre hole



CCI - SC 1022101
pre-slit red PTFE/white Silicone septa, 11mm natural snap-top polypropylene cap, 6mm centre hole



CCI - SC 104101
blue PTFE/white Silicone septa, 11mm natural snap-top polypropylene cap, 6mm centre hole



CCI - SC 1044101
pre-slit blue PTFE/white Silicone septa, 11mm natural snap-top polypropylene cap, 6mm centre hole

1ml Shell Vials - Snap type cap PP star burst Centre. [8.2mm x 40mm] & Inserts



CCI - V181
1ml shell vials, 1st hydrol glass, clear, 8mm PE-Plug, soft, without insertion Residual Volume < 40µl



CCI - I300
300µl Micro-Insert, 31*6mm, clear glass, flat bottom; suits for ND9 Vials Residual Volume < 25µl



CCI - I250
250µl Micro-Insert, 31*5mm, clear glass, flat bottom; suits for ND8 Vials Residual Volume < 20µl



CCI - IP 250
250µl Insert with mandrel interior and polymer feet, 29*5.7mm; suits for ND9 Vials Residual Volume < 2µl



CCI - IP 150
150µl Insert with mandrel interior and polymer feet, 29*5mm; suits for ND8 Vials Residual Volume < 2µl

ND 8 - Screw neck short thread for HPLC / GC. [1.5-2ml of 11.6mm x 32mm]

Vials



CCI - V817
1.5ml small opening screw-thread vial, clear



CCI - V827
1.5ml small opening screw-thread vial with write-on spot, clear



CCI - V835
1.5ml small opening screw-thread vial, amber




CCI - V845
1.5ml small opening screw-thread vial with write-on spot, amber


Caps / Septas




CCI - SC 8281
red PTFE/white Silicone septa, 8mm black screw polypropylene cap, 5.5mm centre hole



CCI - SC 82281
pre-slit red PTFE/white Silicone septa, 8mm black screw polypropylene cap, 5.5mm centre hole



CCI - SC 8481
blue PTFE/white Silicone septa, 8mm black screw polypropylene cap, 5.5mm centre hole



CCI - SC 84481
pre-slit blue PTFE/white Silicone septa, 8mm black screw polypropylene cap, 5.5mm centre hole

ND 13 - Screw neck type Multipurpose. [4ml, 15mm x 45mm]

Vials




CCI - V1317
4ml screw-thread vial, clear




CCI - V1335
4ml screw-thread vial, amber


Caps / Septas



CCI - SC 132131
red PTFE/white Silicone septa, 13mm black screw polypropylene cap, 8.5mm centre hole



CCI - SC 134132
nature PTFE/nature Silicone septa, 13mm black screw polypropylene cap, closed-top



CCI - SC 134131
nature PTFE/nature Silicone septa, 13mm black screw polypropylene cap, 8.5mm centre hole



CCI - C 132
black screw polypropylene cap, closed-top Ø13mm

ND 20 - Head Space Crimp type for HS auto sampler GC [10-20ml, 22.5mm x 75.5mm]

Crimp Type

Vials



CCI - VH1017
10ml crimp headspace clear vial, flat bottom




CCI - VH2013
20ml crimp headspace clear vial, flat bottom


Caps / Septas




CCI - S 203
white PTFE/blue Silicone septa



CCI - SC 201201
white PTFE/white Silicone septa, 20mm crimp-top aluminum cap, 10mm centre hole



CCI - C 201
white PTFE/white Silicone septa



CCI - SC 203201
white PTFE/blue Silicone septa, 20mm crimp-top aluminum cap, 10mm centre hole

ND 18 - Head Space Screw type for magnetic GC. [10-20ml, 22.5mm x 75mm]

Screw Type

Vials




CCI - VA101
10ml clear precision screw-thread vial, round bottom



CCI - VA201
20ml clear precision screw-thread vial, round bottom

Caps / Septas



CCI - SA 001
blue PTFE/white Silicone septa



CCI - SACA 001
blue PTFE/white Silicone septa, 18mm magnetic precision screw-thread metal cap, 8mm centre hole



CCI - CA 001
magnetic precision screw-thread metal cap, 8mm centre hole

ND 22 - Screw neck type for Storage purpose only. [20ml, 40ml, 60 ml] .

Vials



CCI - V2017
20ml screw-thread storage vial, clear



CCI - V2027
20ml screw-thread storage vial, amber



CCI - V4017
40ml screw-thread storage vial, clear



CCI - V4027
40ml screw-thread storage vial, amber



CCI - V6017
60ml screw-thread storage vial, clear


Caps / Septas



CCI - C 222
white screw polypropylene cap, closed-top



CCI - SC 221221
nature PTFE/nature Silicone septa 22*3mm, 22mm white screw polypropylene cap, 15mm centre hole



CCI - C 224
black screw polypropylene cap, closed-top



CCI - SC 221223
nature PTFE/nature Silicone septa 22*3mm, 22mm black screw polypropylene cap, 15mm centre hole

Chromachemie Laboratory Private Limited

No. 103, Model Export Bhavan, 488B, 14th Cross Peenya Industrial Area, Bangalore - 560 058

◆ Tel. : 080 6701 2086 ◆ Email : vialscci@chromachemie.co.in